

# Torque guide

## Recommended tightening torque

Type of product		Torque – Ncm		
Cover Screw		Manual*	Manual*	Manual*
Healing Abutment Healing Abutment Uni ProHeal Cap		Manual*	Manual*	Manual*
TempDesign™ Temporary Abutment		–	15	15
20°/45° Cresco™ Insert for Astra Tech Implant System™ 20°/45° UniAbutment		–	15	15
Bridge Screws Cresco™ Bridge Screw		–	15	15
Atlantis™ for Astra Tech Implant System™ ZirDesign™ TiDesign™ CastDesign™ Cast-to Abutment Angled Abutment		15**	20	25
Direct Abutment Ball Abutment Locator™ Abutment		–	25	25

\*Only light finger force (5–10 Ncm) using a manual screwdriver. Do not use a Ratchet or Torque Wrench

\*\*Note: Only available for TiDesign™